# Cole-Parmer®

Cole-Parmer® Handheld Homogenizer HO-200

- Ergonomic design
- Variable motor speed 5,000 to 35,000rpm
- Processing range from 0.03ml-100ml
- Can be used with robust polycarbonate disposable probes or stainless steel ones







### **Technical Specification**

Specification	HO-200		
Motor	125 watt		
Speed control	5,000 to 35,000rpm		
Processing range	0.03ml to 100ml		
Sound level	<72 db		
Net weight, kg	0.5		
Dimensions, mm (h x Ø)	160 x 55		
Electrical supply	220-240V, 50Hz		
IP Rating	30		

## Cole-Parmer® Handheld Homogenizer HO-200

The Cole-Parmer® HO-200 is a powerful, variable speed, homogenizer. The high-speed, high-torque motor makes the unit ideal for most homogenizing applications. During operation the rotor shaft, which is directly coupled to the motor, spins at up to 35,000rpm. The tube assembly that fits around the shaft remains static and as the rotor shaft spins, within the tube assembly, it creates a pumping action drawing the sample into the open end of the tube, or probe, and forcing it through the windows at the bottom of the tube assembly. This shearing action reduces the particle size of the sample. In addition the speed differential between the rapidly moving rotor, and the static tube assembly creates cavitational force, further dismantling the sample to additionally reduce particle size. The HO-200 is lightweight making it comfortable to use, and most samples can be processed within thirty seconds.

For longer processing the unit can be stand mounted using the included post clamp assembly. A retort stand is available separately. The HO-200 can accommodate a selection of stainless steel probes for larger and tougher samples and also our polycarbonate disposable probes. Disposable probes are suitable for hard tissues and ideal where cross contamination is of concern. If necessary disposable probes can be easily dismantled and autoclaved up to seven times. Through careful selection of probe type, samples of 0.03ml to 100ml can be processed. All Cole-Parmer® homogenizers are supplied with a tool kit for dismantling the rotor probes for easy cleaning

### **Ordering Information**

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Handheld homogenizer, UK plug	04720-01	HO-200	HO-200	SHM1/UK
Handheld homogenizer, European plug	04720-02	HO-200	HO-200-EU	SHM1/EURO

### coleparmer.com

USA: +1.800.323.4340 +1.847.549.7600 Canada: +1.800.363.5900 China: +86.21.5109.9909 France: +33 (0) 1486 37800 Germany: +49 (0) 9377 92030 India: +1.800.266.1244 Italy: +39 (0) 1313 89513 UK: +44 (0) 1480 272279

All other countries: +1.847.549.7600



an antylia scientific company